

Division of Signal Transduction Therapy

Standard Operating Procedure

Preparation of active CAK1 [2 – 368]

<u>Enzyme description:-</u>	CAK1 [2 – 368]
<u>Clone Number:-</u>	DU 1089
<u>Source:-</u>	Recombinant
<u>Expression system:-</u>	<i>E.coli</i>
<u>Tag:-</u>	N-terminal GST and C-terminal His(6) tag
<u>Purification method:-</u>	GSH Sepharose
<u>Expression level:-</u>	<100 µg/L
<u>Calculated molecular mass:-</u>	69, 801 daltons
<u>Purity:-</u>	70 %
<u>Activation protocol:-</u>	Constitutively active
<u>Enzyme storage buffer:-</u>	
	50 mM Tris-HCl pH 7.5, 270 mM Sucrose, 150 mM NaCl, 0.1 mM EGTA, 0.1 % 2-mercaptoethanol, 0.02 % Brij-35, 0.2 mM PMSF, 1 mM Benzamidine.
<u>Storage temperature:-</u>	-70 °C

Assay:-

Two step assay in which CAK1 activates CDK2 [DU 1043] / cyclin A2 [DU 1064]. The activity of CDK2 / cyclin A2 is then assayed against Histone H1 as a substrate (final concentration of 1 mg/ml) in the standard filter binding assay.

Assay buffer:-

50 mM Tris-HCl pH 7.5, 0.1 mM EGTA, 0.1 % 2-mercaptoethanol, 0.1 mM sodium vanadate, 10 mM magnesium acetate

Division of Signal Transduction Therapy

Clone Data Sheet

CAK1 [3 – 368]

<u>Protein</u>	CAK1 [3 – 368]
<u>Clone number</u>	DU 1089
<u>Species</u>	Yeast
<u>Accession number</u>	U60192
<u>Tags</u>	N-terminal GST + C-terminal His(6)
<u>Bacterially expressed protein</u>	<p>MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFELG LEFPNLPYYIDGDVKLTQSMAIIRYIADKHNLGGCPKERAEISMLEGA VLDIRYGVSR IAYS KDFETLKVDFLSKLPEMLKMFEDRLCHKTYLNGDH VTHPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIPOIDKYLKSS KYIAWPLQGWQATFGGGDHPKSDLEVLFGQPLGSKLDSIDITHCQLVK STR TARIYRS DTYAIKCLALDFDIPPHNAKFEVSI LNKLGNKCKHILPL LESKATDNDLLLLFPFEEMNLYEFMQMHYKRDRRKNPYYDLLNPSIP IVADPPVQKYTNQLDVNRYLSFFRQMV EGI AFLHENKIIHRDIKPQNI MLTNTSTVSPKLYIIDFGISYDMANNSQTS AEPMSKVTDISTGIYKA PEVLFGVKCYDGGVDVWSLLIIISQWFQRETSRMGHVPAMIDDGSDDMN SDGSDFRLICSI FEKLGIPSIQKWE EVAQHGSVDAFVGMFGADGDGKYV LDQEKDVQISIVERNMPRLDEIADV KVKQKF INCI LGMV SFSPNERWSC QRILQELEKPGSHHHHHH</p>
<u>Native sequence</u>	<p>Amino acids K2 – P368 (end) of yeast CAK1. Residue K232 of fusion protein is equivalent to K2 of the native enzyme. The GST tag is located at residues 1 - 220 and the His(6) tag is located at residues 601 – 606. The following sequence is present after the CAK1 sequence and before the His(6) tag, GS, residues 599 and 600.</p>
<u>Protease cleavage</u>	PreScission (<u>LEVLFQGPL</u>) at residues 221 – 229
<u>Cloning sites</u>	<i>Bam</i> H1 and <i>Eco</i> R1 site of pGEX6P-1

Division of Signal Transduction Therapy

**Nucleotide
sequence of
insert**

GGATCCAAACTGGATAGTATAGACATTACACACTGTCAGTTG
GTCAAATCTACTAGAAGCTAGGATTTATAGGTCGGATACA
TATGCCATTAAATGTCTAGCACTAGATTTTCGATATCCCGCCA
CATAACGCCAAATTTCGAAGTATCGATATTAAACAAACTGGGC
AACAAATGTAAGCACATCTTACCTCTTCTAGAGTCTAAGGCT
ACCGATAATAATGACCTATTGTTGTTGTTTCCCTTTGAAGAG
ATGAACCTTTATGAGTTCATGCAAATGCACTATAAAAGAGAT
AGAAGAAAAAAAAATCCCTATTACGATTTGCTAAATCCCAGT
ATCCCAATTGTTGCGGACCCCCCGTTCAGAAATATACTAAT
CAATTGGACGTCAATCGGTATTCTTTGTCCTTTTTCCGGCAA
ATGGTTGAAGGGATTGCATTCTTACATGAGAACAAGATCATT
CACCGCGACATCAAACCGCAAATATCATGCTAACAAACAAT
ACCAGCACCGTATCCCCAAAGTTGTACATAATTGATTTTGGC
ATCTCTTATGACATGGCAAATAACTCACAAACAAGTGCGGAA
CCCATGGATAGCAAGGTGACGGATATAAGCACAGGAATTTAC
AAGCCCCAGAAGTGCTTTTTGGAGTAAAATGCTATGATGGT
GGCGTGGACGTGTGGTCGTTGTTGATAATTATTTCTCAGTGG
TTCCAGAGAGAAACAAGCCGTATGGGGCACGTTCCGGCCATG
ATTGATGACGGCAGCGACGACATGAACTCAGATGGAAGCGAT
TTCAGACTGATTTGCTCAATATTCGAAAAGTTGGGCATACCG
TCCATTCAGAAATGGGAAGAGGTTGCGCAACACGGCTCGGTT
GATGCATTTGTTGGTATGTTTGGTGCAGATGGCGATGGCAAG
TATGTACTGGACCAGGAGAAAGATGTACAGATTAGCATTGTT
GAGAGGAATATGCCCTCGACTGGACGAGATTGCGGATGTCAAA
GTCAAGCAGAAGTTCATTAATTGTATCCTGGGGATGGTTTCA
TTTTACCAAACGAAAGATGGAGCTGTCAAAGAATCTTGCAA
GAATTAGAAAAGCCAGGTTACACCACCACCACCACCCTaa
gaattc