

Division of Signal Tranduction Therapy

Standard Operating Procedure

Preparation of active CK1 epsilon [1 – 416]

Enzyme description:- CK1 epsilon [1 – 416]

Clone number:- DU 5127

Source:- Recombinant

Expression system:- *E.coli*

Tag:- N-terminal GST

Purification method:- GSH Sepharose

Calculated molecular mass:-

Monoisotopic 75, 142.42 daltons

Average Mass 75, 190.42 daltons

[cysteines reduced, methionines have not been oxidised

Theoretical pI:- 9.10

Purity:- 85 %

Activation protocol:- Constitutively active

Enzyme storage buffer:-

50 mM Tris-HCl pH 7.5, 270 mM Sucrose, 150 mM NaCl, 0.1 mM EGTA,
0.1 % 2-mercaptoethanol, 0.02 % Brij-35, 0.2 mM PMSF, 1 mM Benzamidine

Storage temperature:- -70 °C

Assay buffer:-

50 mM Tris-HCl pH 7.5, 0.1 % 2-mercaptoethanol, 0.1 mM EGTA, 10 mM MgAc

Substrate:-

KRRRAL S*VASLPGL (where S* is phospho Serine) Final concentration: 300 μM

Division of Signal Tranduction Therapy

Clone Data Sheet

CK1 epsilon [1 – 416]

<u>Protein</u>	CK1 epsilon [1 - 416]
<u>Clone number</u>	DU 5127
<u>Species</u>	Human
<u>Accession number</u>	NM_152221
<u>Tags</u>	N-terminal GST
<u>Bacterially expressed protein</u>	MSPILGYWIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKK FELGLEFPNLPYYIDGDVKLTQSMAIIRYIADKHNLGGCPKERA EISMLEGAVLDIIRYGVSRAYSKDFETLKVDFLSKLPEMLKMFED RLCHKTYLNGDHVTPDFMLYDALDVVLYMDPMCLDAFPKLVCFK KRIEAIPQIDKYLKSSKYIAWPLQGWQATFGGGDHPPKS <u>DLEVLF</u> <u>QGPLGSPEIPGSTRAAAMELRVGNKYRLGRKIGSGSGFGDIYLGAN</u> <u>IASGEEVAIKLECVTKHPQLHIESKFYKMMQGGVGIPSINKCGA</u> EGDYNVMVMELLGPSLEDLFNFCSRKFSLKTVLLLADQMISRIEY IHSKNFIHRDVKPDNFLMGLGKKGNLVYIIDFGLAKKYRDARTHQ HIPYRENKNLTGTARYASINTHLGIEQSRRDDLES LGYVLMYFNL GSLPWQGLKAATKRQKYERISEKKMSTPIEVLCKGYPSEFSTYLN FCRSLRFDDKPDYSYLRLFRNLFHROGFSYDYVFDWNMLKFGAA RNPEDVDRERREHEREERMGQLRGSATRALPPGPPTGATANRLRS AAEPVASTPASRIQPAGNTSPRAISRVDRERKVSMRLHRGAPANV SSSDLTGROEVSRIPASQTSVPFDHLGK
<u>Native sequence</u>	Amino acids M1 – K416 of human CK1 epsilon. Residue M243 of the fusion protein is equilvalent to M1 of the native enzyme. The GST tag is located at residues 1 - 220.
<u>Protease cleavage</u>	PreScission (<u>LEVLFOGP</u>) residues 221 - 228
<u>Cloning sites</u>	<i>NotI</i> sites of pGEX-6P-2

Division of Signal Tranduction Therapy

<u>Nucleotide sequence of insert</u>	gcggccgcgATGGAGCTACGTGTGGGAACAAGTACCG CCTGGGACGGAAGATCGGGAGCGGGTCTTCGGAGATA TCTACCTGGGTGCCAACATCGCCTCTGGTGAGGAAGTC GCCATCAAGCTGGAGTGTGAAGACAAAGCACCCCCA GCTGCACATCGAGAGCAAGTTCTACAAGATGATGCAGG GTGGCGTGGGATCCCGTCCATCAAGTGGTGCAGGAGCT GAGGGCGACTACAACGTGATGGTCATGGAGCTGCTGGG GCCTAGCCTCGAGGACCTGTTCAACTCTGTTCCCGCA AATTCAAGCCTCAAGACGGTGTGCTCTTGGCCGACCAG ATGATCAGCCGCATCGAGTATATCCACTCCAAGAACCTT CATCCACCGGGACGTCAAGCCGACAACCTCCTCATGG GGCTGGGAAAGAAGGGCAACCTGGTCTACATCATCGAC TTCGGCCTGGCCAAGAAGTACCGGGACGCCGCACCCA CCAGCACATTCCCTACCGGGAAAACAAGAACCTGACCG GCACGGCCGCTACGCTTCCATCAACACGCACCTGGGC ATTGAGCAAAGCCGTCGAGATGACCTGGAGAGCCTGGG CTACGTGCTCATGTACTTCAACCTGGCTCCCTGCCCT GGCAGGGGCTCAAAGCAGCCACCAAGGCCAGAAGTAT GAACGGATCAGCGAGAAGAAGATGTCAACGCCATCGA GGTCCTCTGCAAAGGCTATCCCTCCGAATTCTAACAT ACCTCAACTCTGCCGCTCCCTGCCGTTGACGACAAG CCCGACTACTCTTACCTACGTCAGCTCTCCGCAACCT CTTCCACCGCAGGGCTTCTCCTATGACTACGTCTTG ACTGGAACATGCTGAAATTGGTGCAGCCCGGAATCCC GAGGATGTGGACCGGGAGCGGGAGAACACGAACGCGA GGAGAGGATGGGGCAGCTACGGGGTCCCGCACCCGAG CCCTGCCCTGGCCCACCCACGGGGCCACTGCCAAC CGGCTCCGAGTGCCGCCAGCCCCGTGGCTCCACGCC AGCCTCCCGCATCCAGCCGGCTGCCAATACTCTCCA GAGCGATCTCGCGGGTCGACCGGGAGAGGAAGGTGAGT ATGAGGCTGCACAGGGGTGCGCCGCCAACGTCTCCCTC CTCAGACCTCACTGGCGGCAAGAGGTCTCCGGATCC CAGCCTCACAGACAAGTGTGCCATTGACCATCTCGGG AAAGTGA Agcggccgc
--------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------