

Division of Signal Transduction Therapy

Standard Operating Procedure

Preparation of active CLK4 [128 - 481]

<u>Enzyme description:-</u>	CLK4 [128 - 481]
<u>Clone number:-</u>	DU 58744
<u>Source:-</u>	Recombinant
<u>Expression system:-</u>	Baculovirus expression vector system
<u>Tag:-</u>	N-terminal GST
<u>Purification method:-</u>	GSH Sepharose

Calculated molecular mass:-

Monoisotopic 68, 534.60 daltons
Average Mass 68, 578.88 daltons
[cysteines reduced, methionines have not been oxidised]

Theoretical pI:- 6.44

Purity:- >80 %

Activation protocol:- Constitutively active

Enzyme storage buffer:-

50 mM Tris-HCl pH 7.5, 270 mM sucrose, 150 mM NaCl, 0.1 mM EGTA, 10 mM DTT, 0.02 % Brij-35, 0.2 mM PMSF, 1 mM Benzamidine.

Storage temperature:- -70 °C

Assay Buffer:-

50 mM Tris-HCl pH 7.5, 0.1mM EGTA, 10 mM DTT, 10 mM MgAc

Substrate:-

RNRYRDVSPFDHSR Final concentration: 300 uM

Division of Signal Transduction Therapy

Clone Data Sheet

CLK4 [128 - 481]

Protein CLK4 [128 - 481]

Clone number DU 58744

Species Human

Accession number NM_020666.2

Tags N-terminal GST

Baculovirus expressed protein

MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNK
KFELGLEFPNLPYYIDGDVKLTQSMAIIRYIADKHNMLGGCPKE
RAEISMLEGAVLDIRYGVSRIAYSKDFETLKVDFLSKLPMLKM
FEDRLCHKTYLNGDHVTHPDFMLYDALDVVLYMDPMCLDAFPKL
VCFKKRIEAIPOIDKYLKSSKYIAWPLQGWQATFGGGDHPPKSD
LEVLFGQPLGSSKSHRRKRSRSIEDDEEGLICQSGDVLRARYE
IVDTLGEGAFGKVVECIDHGMDGMHVAVKIVKNVGRYREARSE
IQVLEHLNSTDPNSVFRVCVQMLEWFDHHGHVCIVFELLGLSTYD
FIKENSFLPFQIDHIRQMAYQICQSINFLHHNKLTHTDLKPENI
LFVKSDYVVKYNSKMKRDERTLKNTDIKVVDFGSATYDDEHHST
LVSTRHYRAPEVILALGWSQPCDVWSIGCILIEYYLGFVTFQTH
DSKEHLAMMERILGPIPOHMIQKTRKRKYFHHNQLDWDEHSSAG
RYVRRRCKPLKEFMLCHDEEHEKFLDLVRRMLEYDPTQRITLDE
ALQHPFFDLLKKK

Native sequence Amino acids S128 – K481 (end residue) of human CLK4.
Residue S232 of the fusion protein is equivalent to S128 of the
native enzyme. The GST tag is located at residues 1 - 220.

Protease cleavage PreScission (LEVLFQGP) residues 221 - 228

Cloning sites *Bam*H1 and *Not*I sites of pFastBac Dual.

Division of Signal Transduction Therapy

**Nucleotide
sequence of insert**

ggatccTCGAAGAGCCACCGAAGGAAAAGATCCAGGAGTATAGAG
GATGATGAGGAGGGTCACCTGATCTGTCAAAGTGGAGACGTTCTA
AGAGCAAGATATGAAATCGTGGACACTTTGGGTGAAGGAGCCTTT
GGCAAAGTTGTAGAGTGCATTGATCATGGCATGGATGGCATGCAT
GTAGCAGTGAAAATCGTAAAAAATGTAGGCCGTTACCGTGAAGCA
GCTCGTTCAGAAATCCAAGTATTAGAGCACTTAAATAGTACTGAT
CCCAATAGTGTCTTCCGATGTGTCCAGATGCTAGAATGGTTTGAT
CATCATGGTCATGTTTGTATTGTGTTTGAACACTGGGACTTAGT
ACTTACGATTTTCATTAAGAAAACAGCTTTCTGCCATTTCAAATT
GACCACATCAGGCAGATGGCGTATCAGATCTGCCAGTCAATAAAT
TTTTTACATCATAATAAATTAACCCATACAGATCTGAAGCCTGAA
AATATTTTGTGGTGAAGTCTGACTATGTAGTCAAATATAATTCT
AAAATGAAACGTGATGAACGCACACTGAAAACACAGATATCAAA
GTTGTTGACTTTGGAAGTGCAACGTATGATGATGAACATCACAGT
ACTTTGGTGTCTACCCGGCACTACAGAGCTCCCGAGGGTCATTTTG
GCTTTAGGTTGGTCTCAGCCTTGTGATGTTTGGAGCATAGGTTGC
ATTCTTATTGAATATTACCTTGGTTTCACAGTCTTTCAGACTCAT
GATAGTAAAGAGCACCTGGCAATGATGGAACGAATATTAGGACCC
ATACCACAACACATGATTCAGAAAACAAGAAAACGCAAGTATTTT
CACCATAACCAGCTAGATTGGGATGAACACAGTTCTGCTGGTGA
TATGTTAGGAGACGCTGCAAACCGTTGAAGGAATTTATGCTTTGT
CATGATGAAGAACATGAGAACTGTTTGGACCTGGTTCGAAGAATG
TTAGAATATGATCCAACTCAAAGAATTACCTTGGATGAAGCATTG
CAGCATCCTTTCTTTGACTTATTA AAAAAGAAAAtgagcggccgc