

Division of Signal Transduction Therapy

Standard Operating Procedure

Preparation of active PKB gamma [117 - 479]

<u>Enzyme description:-</u>	PKB gamma [117 – 479]
<u>Clone number:-</u>	DU 62934
<u>Source:-</u>	Recombinant
<u>Expression system:-</u>	Baculovirus expression vector system
<u>Tag:-</u>	N-terminal His(6)
<u>Purification method:-</u>	Cobalt Agarose

Calculated molecular mass:-

Monoisotopic 45, 355.23 daltons
Average Mass 45, 384.48 daltons
[cysteines reduced, methionines have not been oxidised]

Theoretical pI:- 5.76

Purity:- >80 %

Activation protocol:- Constitutively active

Enzyme storage buffer:-

50 mM Tris-HCl pH 7.5, 150 mM NaCl, 270 mM sucrose, 0.1 mM EGTA,
0.1 % 2-mercaptoethanol, 0.02 % Brij-35, 1 mM benzamidine, 0.2 mM PMSF

Storage temperature:- -70 °C

Assay buffer:-

50 mM Tris-HCl pH 7.5, 0.1 mM EGTA, 10 mM DTT, 10 mM Magnesium acetate

Substrate:-

GRPRTSSFAEGKK Final concentration: 300 uM

Division of Signal Transduction Therapy

Clone Data Sheet

PKB gamma [117 - 479]

<u>Protein</u>	PKB gamma [117 - 479]
<u>Clone number</u>	DU 62934
<u>Species</u>	Human
<u>Accession number</u>	NM_005465.7
<u>Tags</u>	N-terminal His6
<u>Baculovirus expressed protein</u>	MSYYHHHHHDYDIPTTENLYFQGAMGRSMNCSPTSQIDNIGEEEMDA STTHHKRKTMNDFDYLLKLLGKGTFGKVILVREKASGKYAMKILKKEV IIAKDEVAHTLTESRVLKNTRHPFLTSLKYSFQTKDRLCFVMEYVNGG ELFFHLSRERVFSEDRTRFYGAIEIVSALDYLHSGKIVYRDLKLENLML DKDGHIKITDFGLCKEIGITDAATMKTFCGTPEYLAPEVLEDNDYGRAV DWWGLGVVMYEMMCGRLPFYNQDHEKLFELILMEDIKFPRTLSSDAKS LLSGLLIKDPNKRLGGGPDDAKEIMRHSFFSGVNWQDVYDKKLVPPFK PQVTSETDTRYFDEEFTAQTITITPPEKYDEDGMDCMDNERRPHFPQF SYSASGRE
<u>Native sequence</u>	Amino acids M117 – E479 (end) of human PKB gamma. Residue M30 of the fusion protein is equivalent to M117 of the native enzyme. The His6 tag is located at residues 5 – 10.
<u>Protease cleavage</u>	rTEV (<u>ENLYFQG</u>) residues 18 – 24
<u>Cloning sites</u>	<i>Bgl</i> III <i>Eco</i> R1 into <i>Bam</i> H1 <i>Eco</i> R1 sites into pFastBac HTb

Division of Signal Transduction Therapy

**Nucleotide
sequence of insert**

ATGAATTGTAGTCCAACCTTCACAAATTGATAATATAGGAGAGGAAGAG
ATGGATGCCTCTACAACCCATCATAAAAGAAAGACAATGAATGATTTT
GACTATTTGAACTACTAGGTAAAGGCACTTTTGGGAAAGTTATTTTG
GTTTCGAGAGAAGGCAAGTGGAAAATACTATGCTATGAAGATTCTGAAG
AAAGAAGTCATTATTGCAAAGGATGAAGTGGCACACACTCTAACTGAA
AGCAGAGTATTAAGAACACTAGACATCCCTTTTTTAACATCCTTGAAA
TATTCCTTCCAGACAAAAGACCGTTTTGTGTTTTGTGATGGAATATGTT
AATGGGGCGAGCTGTTTTTCCATTTGTCGAGAGAGCGGGTGTCTCT
GAGGACCGCACACGTTTTCTATGGTGCAGAAATTGTCTCTGCCTTGAC
TATCTACATTCGGAAAGATTGTGTACCGTGATCTCAAGTTGGAGAAT
CTAATGCTGGACAAAGATGGCCACATAAAAATTACAGATTTTGGACTT
TGCAAAGAAGGGATCACAGATGCAGCCACCATGAAGACATTCTGTGGC
ACTCCAGAATATCTGGCACCAGAGGTGTTAGAAGATAATGACTATGGC
CGAGCAGTAGACTGGTGGGGCCTAGGGTTGTCATGTATGAAATGATG
TGTGGGAGGTTACCTTTCTACAACCAGGACCATGAGAACTTTTTGAA
TTAATATTAATGGAAGACATTAATTTCTCGAACACTCTCTTCAGAT
GCAAATCATTGCTTTTCAGGGCTTTGATAAAGGATCCAATAAACGC
CTTGGTGGAGGACCAGATGATGCAAAGAAATTATGAGACACAGTTTC
TTCTCTGGAGTAACTGGCAAGATGTATATGATAAAAAGCTTGTACCT
CCTTTTAAACCTCAAGTAACATCTGAGACAGATACTAGATATTTTGAT
GAAGAATTTACAGCTCAGACTATTACAATAACACCACCTGAAAAATAT
GATGAGGATGGTATGGACTGCATGGACAATGAGAGGCGCCGCATTTTC
CCTCAATTTTCTACTCTGCAAGTGGACGAGAAaagaattc