

Division of Signal Transduction Therapy

Standard Operating Procedure

Preparation of active DAPK2 [1 – 370]

Enzyme description:- DAPK2 [1 - 370]

Clone number:- DU 31115

Source:- Recombinant

Expression system:- *E.coli*

Tag:- N-terminal GST

Purification method:- GSH Sepharose

Calculated molecular mass:-

Monoisotopic 69,677.97 daltons

Average Mass 69,722.26 daltons

[cysteines reduced, methionines have not been oxidised]

Theoretical pI:- 6.16

Purity:- >80 %

Activation protocol:- Constitutively active

Enzyme storage buffer:-

50 mM Tris-HCl pH 7.5, 270 mM Sucrose, 150 mM NaCl, 0.1 mM EGTA,
0.1 % 2-mercaptoethanol, 0.02 % Brij-35, 1 mM benzamidine, 0.2 mM PMSF

Storage temperature:- -70 °C

Assay buffer:-

50 mM Tris-HCl pH 7.5, 0.1 mM EGTA, 0.1 % 2-mercaptoethanol, 10 mM MgAc,

Substrate:-

KKLNRTL SFAEPG

Final concentration: 300 μ M

Division of Signal Transduction Therapy

Clone Data Sheet

DAPK2 [1 – 370]

Protein DAPK2 [1 - 370]

Clone number DU 31115

Species Human

Accession number NM_014326.3

Tags N-terminal GST

Bacterially expressed protein

MSPILGYWKIKGLVQPTRLLLEYLEEKYEHLIERDEGDKWRNKKFELG
LEFPNLPYYIDGDVKLTQSMAIIRYIADKHNMLGGCPKERAETSMLEGA
VLDIRYGVSR IAYS KDFETLKVDFLSKLP EMLKMFEDRLCHKTYLNGDH
VTHPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIPOIDKYLKSS
KYIAWPLQGWQATFGGGDHPPKSDLEVLFGQPLGSMFQASMRSPNMEPF
KQQKVEDFYDIGEELGSGQFAIVKKCREKSTGLE YAAFKIKKRQSRASR
RGVSREEIEREVSILRQVLHNVITLHDVYENRTDVVLILELVSGGELF
DFLAQKESLSEEEATSF IKQILDGVNYLHTKKIAHF DLKPENIMLLDKN
IPIPHIKLIDFGLAHEIEDGVEFKNIFGTPEFVAPEIVNYEPLGLEADM
WSIGVITYILLSGASPFLGDTKQETLANITAVSYDFDEEFFSQTSELAK
DFIRKLLVKETRKR LTIQEALRHPWITPVDNQQAMVRRESVVNLENFRK
QYVRRRWKLSFSIVSLCNHLTRSLMKKVHLRPDEDLRNCESDTEEDIAR
RKALHPRRRSSTS

Native sequence Amino acids M1 – S370 (end) of human DAPK2.
Residue M232 of the fusion protein is equivalent to M1 of the native enzyme. The GST tag is located at residues 1 – 220.

Protease cleavage PreScission (LEVLFQGP) residues 221 - 229

Cloning sites *Bam*H1 and *Not*I sites of pGEX 6P-1

Division of Signal Transduction Therapy

Nucleotide
Sequence of insert

ggatccATGTTCCAGGCCTCAATGAGGAGTCCAAACATGGAGCCATTCA
AGCAGCAGAAGGTGGAGGACTTTTATGACATCGGAGAGGAGCTGGGGAG
TGGCCAGTTTGCCATCGTGAAGAAGTGCCGGGAGAAGAGCACGGGGCTT
GAGTATGCAGCCAAGTTCATCAAGAAGCGGCAGAGCCGGGCGAGCCGGC
GCGGTGTGAGCCGGGAGGAGATCGAGCGGGAGGTGAGCATCCTGCGGCA
GGTGCTGCACCACAATGTCATCACGCTGCACGACGTCTATGAGAACCGC
ACCGACGTGGTGCTCATCCTTGAGCTAGTGTCTGGAGGAGAGCTCTTCG
ATTCCTGGCCCAGAAGGAGTCACTGAGTGAGGAGGAGGCCACCAGCTT
CATTAAGCAGATCCTGGATGGGGTGAACCTCACACAAAGAAAATT
GCTCACTTTGATCTCAAGCCAGAAAACATTATGTTGTTAGACAAGAATA
TTCCCATTCACACATCAAGCTGATTGACTTTGGTCTGGCTCACGAAAT
AGAAGATGGAGTTGAATTTAAGAATATTTTTGGGACGCCGGAATTTGTT
GCTCCAGAAATTGTGAACTACGAGCCCCTGGGTCTGGAGGCTGACATGT
GGAGCATAGGCGTCATCACCTACATCCTCTTAAGTGGAGCATCCCCTT
CCTGGGAGACACGAAGCAGGAAACACTGGCAAATATCACAGCAGTGAGT
TACGACTTTGATGAGGAATTCCTTCAGCCAGACGAGCGAGCTGGCCAAG
ACTTTATTCGGAAGCTTCTGGTTAAAGAGACCCGGAACGGCTCACAAT
CCAAGAGGCTCTCAGACACCCCTGGATCACGCCGGTGGACAACCAGCAA
GCCATGGTGCGCAGGGAGTCTGTGGTCAATCTGGAGAACTTCAGGAAGC
AGTATGTCCGCAGGCGGTGGAAGCTTTCCTTCAGCATCGTGTCCCTGTG
CAACCACCTCACCCGCTCGCTGATGAAGAAGGTGCACCTGAGGCCGGAT
GAGGACCTGAGGAACTGTGAGAGTGACACTGAGGAGGACATCGCCAGGA
GGAAAGCCCTCCACCCACGGAGGAGGAGCAGCACCTCCtaagcggccgc