

MRC PPU Reagents and Services

Standard Operating Procedure

Preparation of active TNK2 [1 - 389]

<u>Enzyme description:-</u>	TNK2 [1 – 389]
<u>Clone number:-</u>	DU 60935
<u>Source:-</u>	Recombinant
<u>Expression system:-</u>	Baculovirus expression vector system
<u>Tag:-</u>	N-terminal GST
<u>Purification method:-</u>	GSH Sepharose

Calculated molecular mass:-

Monoisotopic	71, 334.53 daltons
Average Mass	71, 380.55 daltons

[cysteines reduced, methionines have not been oxidised]

Theoretical pI:- 7.03

Purity:- >75 %

Activation protocol:- Constitutively Active

Enzyme storage buffer:-

50 mM Tris-HCl pH 7.5, 270 mM Sucrose, 150 mM NaCl, 0.1 mM EGTA, 0.1 % 2-mercaptoethanol, 0.02 % Brij-35, 0.2 mM PMSF, 1 mM Benzamidine.

Storage temperature:- -70 deg C

Assay buffer:-

50 mM Tris-HCl pH 7.5, 10 mM DTT, 10 mM magnesium acetate, 0. 1mM EGTA

Substrate:-

EAIYAAPFAKKK Final concentration: 300 µM

MRC PPU Reagents and Services

Clone Data Sheet

TNK2 [1 – 389]

Protein TNK2 [1 – 389]

Clone number DU 60935

Species Human

Accession number NM_005781

Tags N-terminal GST

Baculovirus expressed protein
MSPILGYWKIKGLVQPTRLLEYLEEKYEEHLYERDEGDKWRNKKFEL
GLEFPNLPHYIDGDVKLTSMAIIRYIADKHNMLGGCPKERAIEISMLE
GAVLDIRYGVSRIAYSKDFETLKVDFLSKLPPEMLKMFEDRLCHKTYLN
GDHVTHPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIPOIDKY
LKSSKYIAWPLQGWQATFGGGDHPKSDLEVLFOGPLGSMQPEEGTGW
LLELLSEVQLQOYFLRLRDDLNVTRLSHFYVKNEDLEKIGMGRPGQR
RLWEAVKRRKALCKRKSWSKVFSGKRLEAEFPPHHSQSTFRKTSPPAP
GGPAGEGPLQSLTCLIGEKDLRLLLEKLGDSFGVRRGEWDAPSGKTV
SVAVKCLKPDVLSQPEAMDDFIREVNAMHSLDHRNLIIRLYGVVLTTPM
KMVELAPLGSLLDRLRKHQGHFLLGTLSRYAVQVAEGMGYLESKRFI
HRDLAARNLLLATRDLVKIGDFGLMRALPQNDHYVMQEHRKVPFWC
APESLKTRTFSHASDTWMFGVTLWEMFTYGOEPWIGLNGSQILHKIDK
EGERLPRPEDCPQDIYNVMVQCWAHKPEDRPTFVALRDFLLEAQ

Native sequence Amino acids M1 – Q389 (end residue is R1038) of human TNK2.
Residue M232 of the fusion protein is equivalent to M1 of the native
enzyme. The GST tag is located at residues 1 – 220.

Protease cleavage PreScission (LEVLFQGP) residues 221 - 228

Cloning sites *Bam*H1 - *Not*I sites of pFastBac GST

MRC PPU Reagents and Services

Nucleotide
sequence of
insert

ggatccATGCAGCCAGAGGAGGGCACAGGCTGGCTGCTGGAGCTGCTG
TCCGAGGTGCAGCTGCAACAGTACTTCCTGCGGCTCCGAGATGACCTC
AACGTCAACCGCCTGTCCCACCTTTGAGTACGTCAAGAATGAGGACCTG
GAGAAGATCGGCATGGGTGCGCCTGGCCAGCGGCGGCTGTGGGAGGCT
GTGAAGAGGAGGAAGGCCCTTGTGCAAACGCAAGTCGTGGATGAGTAAG
GTGTTCACTGGAAAGCGACTGGAGGCTGAGTTCCCACCTCATCACTCT
CAGAGCACCTTCCGGAAGACCTCGCCCCGCCCTGGGGGCCAGCAGGG
GAGGGGCCCTGCAGAGCCTCACCTGCCTCATTGGGGAGAAGGACCTG
CGCCTCCTGGAGAAGCTGGGTGATGGTTCCCTTTGGCGTGGTGCGCAGG
GGCGAGTGGGACGCGCCCTCAGGGAAGACGGTGAGTGTGGCTGTGAAG
TGCCTGAAGCCCGATGTCCTGAGCCAGCCAGAAGCCATGGACGACTTC
ATCCGGGAGGTCAATGCCATGCACTCGCTCGACCACCGAAACCTCATC
CGCCTCTACGGGGTGGTGCTCACGCCGCCATGAAGATGGTGACAGAG
CTGGCACCTCTGGGATCGTTGTTGGACCGGCTACGTAAGCACCAGGGC
CACTTCCTCCTGGGACTCTGAGCCGCTACGCTGTGCAGGTGGCTGAG
GGCATGGGCTACCTGGAGTCCAAGCGCTTTATTCACCGTGACCTGGCT
GCCCCGAATCTGCTGTTGGCTACCCGCGACCTGGTCAAGATCGGGGAC
TTTGGGCTGATGCGAGCACTACCTCAGAATGACGACCATTACGTCATG
CAGGAACATCGCAAGGTGCCCTTCGCCTGGTGTGCCCCGAGAGCCTG
AAGACACGCACCTTCTCCCATGCCAGCGACACCTGGATGTTTCGGGGTG
ACACTGTGGGAAATGTTACCTACGGCCAGGAGCCCTGGATCGGCCTC
AACGGCAGTCAGATCCTGCATAAGATCGACAAGGAGGGGGAGCGGCTG
CCCCGGCCCGAGGACTGTCCCAGGACATCTACAACGTCATGGTCCAG
TGCTGGGCTCACAAAGCCAGAGGACAGACCCACGTTTGTGGCCCTGCGG
GACTTCCTGCTGGAGGCCAGtgagcggccgc