

MRC PPU Reagents and Services

Standard Operating Procedure

Preparation of Envelope Protein [1 – 82] MERS CoV

Enzyme description:- MERS E Protein [1 - 82]

Clone number:- DU 68503

Source:- Recombinant

Expression system:- *E. coli*

Tag:- N-terminal GST

Purification method:- GSH Agarose

Calculated molecular mass:-

Monoisotopic 36, 154.49 daltons

Average Mass 36, 178.40 daltons

[cysteines reduced, methionines have not been oxidised]

Theoretical pI:- 5.92

Purity:- 65 %

Enzyme storage buffer:-

50 mM Tris-HCl pH 7.5, 270 mM Sucrose, 150 mM NaCl, 0.1 mM EGTA,
0.1 % 2-mercaptoethanol, 0.03 % Brij-35

Storage temperature:- -70 °C

MRC PPU Reagents and Services

Clone Data Sheet

Envelope Protein [1 – 82] MERS CoV

<u>Protein</u>	MERS E Protein [1 - 82]
<u>Clone number</u>	DU 68503
<u>Accession number</u>	AGV08472.1
<u>Tags</u>	N-terminal GST
<u>Bacterially expressed protein</u>	<p>MSPILGYWKIKGLVQPTRLLLEYLEEKYEEHLYERDEGDKWRNKK FELGLEFPNLPYYIDGDVKLTQSMAIIRYIADKHNMLGGCPKERA EISMLEGAVLDIRYGVSRIAYSKDFETLKVDFLSKLPPEMLKMFED RLCHKTYLNGDHVTHPDFMLYDALDVVLYMDPMCLDAFPKLVCFK KRIEAIPOIDKYLKSSKYIAWPLQGWQATFGGDHPPKSDLEVL <u>QGPLGSMLPFVQERIGLFI</u>VNFFIFTVVCAITLLVCM<u>AF</u>L<u>TATRL</u> CVQCMTGFNTLLVQPALYLYNTGRSVYVKFQDSKPPLPPDEWV</p>
<u>Native sequence</u>	<p>Amino acids M1 – V82 (end). Residue M232 of the fusion protein is equivalent to M1 of the native enzyme. The GST tag is located at residues 1 – 220.</p>
<u>Protease cleavage</u>	PreScission (<u>LEVLFQGP</u>) residues 221 – 228